

## Specifications

<b>Model</b>	HEM-7361T		
<b>Features</b>	 Intelli Wrap™	 AFIB Indicator	 connect
	 Dual Users	 100 Memory	 Intelli sense
	 Irregular Heartbeat Detection	 Average Of Last 3 Readings	 Hypertension Indicator
	 Cuff Wrapping Guide	 Pulse	 Body Movement Detection
<b>Measurement Method</b>	Oscillometric Method		
<b>Blood Pressure Measurement Range</b>	SYS: 60 to 260 mmHg DIA: 40 to 215 mmHg		
<b>Pulse Measurement Range</b>	40 to 180 beats/min		
<b>Accuracy</b>	Pressure: ±3 mmHg   Pulse: ±5% of reading		
<b>Battery Life</b>	Capacity of new alkaline batteries is approx. 1000 measurements		
<b>Power Source</b>	4 "AA" alkaline batteries 1.5V or AC adapter (optional, INPUT AC100-240V 50-60Hz)		
<b>Console Weight</b>	Approximately 460g without batteries		
<b>Outer Dimensions</b>	Approximately 191 mm (W) × 85mm (H) × 120 mm (L)		
<b>Arm Cuff</b>	Fits arm circumference 22 - 42cm		
<b>Package Content</b>	Monitor, Intelli Wrap™ Cuff, storage case, 4 AA batteries, instruction manual (1) and (2), OMRON connect setup instructions		

<sup>1</sup> Odotayo A, Wong CX, Hsiao AJ, Hopewell S, Altman DG et al. (2016) Atrial fibrillation and risks of cardiovascular disease, renal disease, and death: systematic review and meta-analysis. *Bmj* 354:i4482.

<sup>2</sup> Camm AJ, Lip GY, De Caterina R, Savelieva I, Atar D, Hohnloser SH et al.: 2012 focused update of the ESC Guidelines for the management of atrial fibrillation: an update of the 2010 ESC Guidelines for the management of atrial fibrillation--developed with the special contribution of the European Heart Rhythm Association. *Europace* 2012; 14:1385-1413.

<sup>3</sup> Hart RG, Benavente O, McBride R, Pearce LA: Antithrombotic therapy to prevent stroke in patients with atrial fibrillation: a meta-analysis. *Ann Intern Med* 1999; 131:492-501.

<sup>4</sup> Ruff CT, Giugliano RP, Braunwald E, Hoffman EB, Deenadayalu N, Ezekowitz MD et al.: Comparison of the efficacy and safety of new oral anticoagulants with warfarin in patients with atrial fibrillation: a meta-analysis of randomised trials. *Lancet* 2014; 383:955-962.

<sup>5</sup> Morillo, Carlos A et al. "Atrial fibrillation: the current epidemic." *Journal of geriatric cardiology : JGC* vol. 14,3 (2017): 195-203. doi:10.11909/j.issn.1671-5411.2017.03.011.

<sup>6</sup> NHS, Atrial Fibrillation, 2018. Available at <https://www.nhs.uk/conditions/atrialfibrillation/complications/>. Last accessed May 2019.

<sup>7</sup> Camm AJ, Savelieva I, Bharati S, et al. Atrial tachycardia, flutter, and fibrillation. In: Saksena S, Camm AJ, eds. *Electrophysiological Disorders of the Heart*. Elsevier Churchill- Livingstone, 2005:283-363.

### Intended Use

This device is a digital monitor intended for use in measuring blood pressure and pulse rate in adult patient population. The device detects the appearance of irregular heartbeats during measurement and gives a warning signal with readings. It is mainly designed for general household use.

The device can detect an irregular pulse suggestive of Atrial Fibrillation (Afib). Please note that the device is not intended to diagnose Afib. A diagnosis of Afib can only be confirmed by Electrocardiogram (ECG). If the Afib symbol appears, consult your physician.

Refer to the bottom panel for manufacturer information.

This OMRON product is produced under the strict quality system of OMRON HEALTHCARE Co., Ltd., Japan.


The core component for OMRON blood pressure monitors, which is the Pressure Sensor, is produced in Japan.

# YOUR BEST CHANCE AT PREVENTING A STROKE

1-in-4 People Are At Risk Of Afib Related Strokes

# JUST GOT BETTER

**AFIB ALERT!**



**OMRON**




\*Based on the results of an independent quantitative survey conducted by MIMS in Jul-Sep2019 on total 600 doctors (including cardiologists, physicians and specialists) in India, Indonesia, Philippines, Malaysia and Vietnam who have recommended BPM products to their patients.

# EARLY DETECTION ENABLES EARLY TREATMENT


Early detection of Atrial Fibrillation (AFib) followed by adequate treatment can reduce the risk of a stroke by 68%<sup>1,2,3</sup>. Omron HEM-7361T with AFib detection is the best first step to detect any abnormalities in your heart rhythm and in your blood pressure, from the comfort of home.

## AFIB BY THE NUMBERS



## STAY ON TOP OF YOUR HEART HEALTH WITH HEM-7361T DUAL CHECK: BP MONITORING + AFIB ALERT

Detect atrial fibrillation while measuring your blood pressure in the comfort of your home, at anytime with the HEM-7361T. Equipped with breakthrough AFib detection technology, among other advanced features, you will instantly be alerted of any AFib abnormalities in heart rhythm, allowing you to consult doctors for early treatments.



## DUAL CHECK FOR DUAL USERS


The OMRON HEM-7361T measures and detects high blood pressure and AFIB at the same time. This allows users to have a comprehensive check for their risk of stroke with one device, in one sitting. Designed with families in mind, the dual user function allows for the storage and transmission of measurement data of 2 separate users to their smart phones.

## SEAMLESS SYNC TO OMRON CONNECT

With the Bluetooth® enabled feature of the OMRON HEM-7361T, you can now seamlessly sync your daily health data to your OMRON Connect profile for a better overview of your health, not just for you, but for your doctor as well via email.

## CLINICALLY VALIDATED

OMRON Blood Pressure Monitors are most trusted by Doctors for their heart health monitoring as OMRON Blood Pressure Monitors have been clinically validated and meet the highest standards of precision and reliability from trusted organisations such as Association for the Advancement of Medical Instrumentation Standard (AAMI) and European Society for Hypertension (ESH).



OMRON



## KEY FEATURES

- AFIB Indicator
- Hypertension Indicator
- Intelli Wrap™
- Dual Users
- Irregular Heartbeat Detection
- OMRON Connect
- Intelli sense
- Dual Check

